

***In Vitro* Generation of Antigen-inexperienced Memory CD8⁺ T Cells by IL-15**

Background

Memory-phenotype CD8⁺ T cells exist even within hosts that have not been exposed to foreign antigen (e.g., in human cord blood).¹ Classically, memory CD8⁺ T cells (T_{MEM}) result from antigen-specific interactions, but immunological memory can also be achieved through antigen-inexperienced memory-phenotype T cells — virtual memory CD8⁺ T cells (T_{VM}). T_{VM} arise in the periphery from naive CD8⁺ T cells and are driven by homeostatic cytokines.^{1,2} They partially resemble classical T_{MEM} with their ability to rapidly mount an immune response upon antigen stimulation.^{1,2,3} Putative human T_{VM} express NK and T_{MEM} markers (EOMES, NKG2A, KIR3DL1) and are phenotypically defined as [CD45RA+CD27-EOMES+NKG2A+panKIR+]⁹. It has been shown that T_{VM} are largely dependent on IL-15 and IL-4 for differentiation, expansion, and survival.^{1,4,5} Ablation of IL-15 production or IL-15 signaling completely inhibits T_{VM} differentiation and survival.⁵ BATF-3 dependent dendritic cells (DC) have been shown to be the predominant producers of IL-15 in the steady state, and BATF-3 deficient mice had a severely impaired T_{VM} population. T_{VM} differentiation in IL-15 deficient mice can be rescued by injecting IL-15/IL-15RA complexes.⁶ Ablation of IL-4 is also known to compromise the T_{VM} population, indicating their dependence on IL-4.^{6,7,8} However, this cytokine dependency has mainly been shown in mice. Whether IL-15 and/or IL-4 directly contributes to T_{VM} differentiation in humans is to be determined.

Hypothesis

Treatment of bulk CD8⁺ naive T cells with IL-15 and/or IL-4 *in vitro* will convert them into [CD45RA+CD27-EOMES+panKIR+] antigen-inexperienced memory CD8⁺ T cells (T_{VM}).

Rationale

It has been shown that CD45RA+CD27-panKIR+ T_{VM} cells accumulate in patients with human immunodeficiency virus (HIV). They are positively correlated with IL-15 and IFN- γ serum concentrations.⁹ However, there is still no clear evidence for cytokine-induced activation of T_{VM} in human infectious diseases. T_{VM} are also inversely correlated with latent HIV viral load and have been demonstrated to play a protective role in influenza, helminth infections, and other chronic infections.^{10,11,12}

A putative T_{VM} equivalent population in humans has been defined. Given that human T_{VM} have a similar responsiveness to antigen stimulation and phenotype to mouse T_{VM}, and that they are positively correlated with IL-15 serum levels in patients with chronic infections, it is plausible that IL-15 mediates the expansion/differentiation of human T_{VM} cells as well.^{1,2,9}

Objectives

1. Establish proportion of Tvm in isolated Tn / bulk CD8⁺
2. Generate Tvm from bulk CD8⁺ T cells
3. Generate Tvm from isolated Tn
4. Assess the function of Tvm following 48hrs CD3/28 stimulation

Methods

Establish proportion of Tvm in PBMCs

1. Using flow cytometry, sort PBMCs with CD3,CD8,Livedead,CD45RA,CD27,CCR7, EOMES, KIR3DL1,KIR2D,NKG2A. Determine the Tvm population percentage, gated on CD3+CD8+CD45RA+CCR7-CD27-EOMES+KIR3DL1+ and/or KIR2D+ and NKG2A+⁵

Generate Tvm from PBMC culture

1. Thaw PBMCs in a 6-well plate overnight (or 24hrs)
2. Coculture PBMCs in round-bottom well with:
a. control (no cytokines)
b. IL-15 (20ng/mL)
c. IL-4 (20ng/mL)
d. IL-15 + IL-4 ((20ng/mL each)
3. Using flow cytometry, sort groups a-d with CD3,CD8,Livedead,CD45RA,CD27,CCR7, EOMES, KIR3DL1,KIR2D,NKG2A. Determine the Tvm population percentage, gated on CD3+CD8+CD45RA+CCR7-CD27-EOMES+KIR3DL1+ and/or KIR2D+ and NKG2A+⁵

Generate Tvm from bulk CD8+ T cells

1. Thaw PBMCs in a 6-well plate overnight (or 24hrs)
2. Isolate CD8+ T cells using CD8 sorting kit.
3. Coculture bulk CD8 T cells with:
a. control (no cytokines)
b. IL-15 (20ng/mL)
c. IL-4 (20ng/mL)
d. IL-15 + IL-4 ((20ng/mL each)
4. Sort cells by flow cytometry according to the panel in (3) above, w/out CD3

Generate Tvm from isolated Tn

1. Thaw PBMCs in a 6-well plate overnight (or 24hrs)
2. Isolate CD8+CD45RA+CCR7+ Tn using Tn magnetic sorting kit.
3. Coculture Tn cells with:
a. control (no cytokines)
b. IL-15 (20ng/mL)
c. IL-4 (20ng/mL)
d. IL-15 + IL-4 ((20ng/mL each)
4. Sort cells by flow cytometry according to the panel in (3) above, w/out CD3

****Putative human Tvm markers are CD8+CD3+CD45RA+CD45RO-CD27-CCR7-EOMES+KIR3DL1+ and/or KIR2D and NKG2A+.^{5,9} KIR3DL1 and KIR2D collectively comprise the panKIR marker, which may be the ideal marker for Tvm Identification, rather than EOMES.¹⁷**

Assess the function of Tvm following 48hrs CD3/28 stimulation

1. Isolate CD8+CD45RA+CCR7+ Tn using Tn magnetic sorting kit.

2. Coculture CD45RA+CCR7+ bulk Tn with:

a. IL-15+IL-4+IL-2

b. IL-15+IL-2

c. IL-4+IL-2

d. Control (IL-2)

*10ng/mL each cytokine for 7 days^{4,16}

3. For groups a-d, stimulate with CD3+CD28 for 48hrs and split into another unstim population:

Treatment, Stim (ts)

as. IL-15+IL-4+IL-2; CD3/28

bs. IL-15+IL-2; CD3/28

cs. IL-4+IL-2; CD3/28

Treatment, Unstim (tu)

au. IL-15+IL-4+IL-2

bu. IL-15+IL-2

cu. IL-4+IL-2

Control, Stim (cs)

ds. IL-2; CD3/28

Control, Unstim (cu)

du. IL-2

4. Using flow cytometry, sort above groups with CD8,Livedead,CD45RA,CD27,CCR7,EOMES,panKIR, NKG2A,granzyme B,perforin,(granulysin).⁹ Compare granzyme B and perforin expression between ts and cs.

References

1. White, J. T., Cross, E. W. & Kedl, R. M. Antigen-inexperienced memory CD8+ T cells: where they come from and why we need them. *Nat Rev Immunol* 17, 391–400 (2017).

2. Hussain, T. & Quinn, K. M. Similar but different: virtual memory CD8 T cells as a memory like cell population. *Immunol Cell Biol* 97, 675–684 (2019).

3. Quinn, K. M. *et al.* Age-Related Decline in Primary CD8+ T Cell Responses Is Associated with the Development of Senescence in Virtual Memory CD8+ T Cells. *Cell Reports* 23, 3512–3524 (2018).

4. Alves, N. L., Hooibrink, B., Arosa, F. A. & van Lier, R. A. W. IL-15 induces antigen-independent expansion and differentiation of human naive CD8⁺ T cells in vitro. *Blood* 102, 2541–2546 (2003).
5. White, J. T. *et al.* Virtual memory T cells develop and mediate bystander protective immunity in an IL-15-dependent manner. *Nat Commun* 7, 11291 (2016).
6. Sosinowski, T. *et al.* CD8 α ⁺ Dendritic Cell *Trans* Presentation of IL-15 to Naive CD8⁺ T Cells Produces Antigen-Inexperienced T Cells in the Periphery with Memory Phenotype and Function. *J.I.* **190**, 1936–1947 (2013).
7. Akue, A. D., Lee, J.-Y. & Jameson, S. C. Derivation and Maintenance of Virtual Memory CD8 T Cells. *J.I.* **188**, 2516–2523 (2012).
8. Kurzweil, V., LaRoche, A. & Oliver, P. M. Increased Peripheral IL-4 Leads to an Expanded Virtual Memory CD8⁺ Population. *J.I.* **192**, 5643–5651 (2014).
9. Jin, J.-H. *et al.* Virtual memory CD8⁺ T cells restrain the viral reservoir in HIV-1-infected patients with antiretroviral therapy through derepressing KIR-mediated inhibition. *Cell Mol Immunol* 17, 1257–1265 (2020).
10. Quinn, K. M. & Hussain, T. Bystanders or real players: virtual memory T cells keep chronic infections in check. *Cell Mol Immunol* 17, 797–798 (2020).
11. Lee, H., Jeong, S. & Shin, E.-C. Significance of bystander T cell activation in microbial infection. *Nat Immunol* 23, 13–22 (2022).
12. Lin, J. S. *et al.* Virtual memory CD8 T cells expanded by helminth infection confer broad protection against bacterial infection. *Mucosal Immunol* 12, 258–264 (2019).
13. Webb, G. M. *et al.* The human IL-15 superagonist ALT-803 directs SIV-specific CD8⁺ T cells into B-cell follicles. *Blood Advances* 2, 76–84 (2018).
14. Alanio, C. *et al.* Bystander hyperactivation of preimmune CD8⁺ T cells in chronic HCV patients. *eLife* 4, e07916 (2015).
15. Kakumu, S. *et al.* Serum levels of IL-10, IL-15 and soluble tumour necrosis factor-alpha (TNF- α) receptors in type C chronic liver disease. *Clinical and Experimental Immunology* **109**, 458–463 (2003).
16. Geginat, J., Lanzavecchia, A. & Sallusto, F. Proliferation and differentiation potential of human CD8⁺ memory T-cell subsets in response to antigen or homeostatic cytokines. *Blood* 101, 4260–4266 (2003).
17. Jacomet, F. *et al.* Evidence for eomesodermin-expressing innate-like CD8⁺ KIR/NKG2A⁺ T cells in human adults and cord blood samples: Innate immunity. *Eur. J. Immunol.* 45, 1926–1933 (2015).